



E-Prescription Recommendation

Internship at Prixa.ai

Dave Danadiva Agusta Peerera (12006002)

dave.peerera@student.sgu.ac.id

INTRODUCING



01

VISION

To enable a promising future for you and the next generation

02

MISSION

Humanize healthcare using data and technology

WEB-APP SERVICES

Symptom Checker 01

Teleconsultation 02


Lab Test 03

Pharmacy Delivery 04

Doctor Appointment 05

For more info visit: prixa.ai

PRODUCT TEAM STRUCTURE



E-Prescription Recommendation




Goals

Develop and evaluate a medicine drug recommendation system that can assist doctors in making decision faster, more precise, and accurate prescriptions for patients with certain diseases.




BACKGROUND

Doctor spent approximately 23 minutes for each teleconsultation session




BACKGROUND

Doctor spent approximately **23 minutes** for each teleconsultation session



TOOLS



Python 3.7 & Python 3.10



Google Colaboratory IDE



Visual Studio Code IDE





Database



Visualization

METHODOLOGY PROJECT STAGES





STAGE 1 METHODOLOGY

DATA GATHERING & CLEANING

- Gather data 01
- Convert data file format 02
- Filter out dummy data 03
- Filter out rubbish data 04

01

Explore Data gathered to gain insight and obtain correlation between features in the dataset

02

Exploratory data analysis will be divided into patients and prescriptions part and medicine part




STAGE 2

METHODOLOGY

**EXPLORATORY
DATA ANALYSIS**

METHODOLOGY


EXPLORATORY DATA ANALYSIS



PART I: PATIENTS & PRESCRIPTIONS

METHODOLOGY

EXPLORATORY DATA ANALYSIS




PART I: PATIENTS & PRESCRIPTIONS

METHODOLOGY

EXPLORATORY DATA ANALYSIS


Patient Distribution by Life Stages



PART I: PATIENTS & PRESCRIPTIONS

METHODOLOGY


EXPLORATORY DATA ANALYSIS



PART I: PATIENTS & PRESCRIPTIONS

METHODOLOGY


EXPLORATORY DATA ANALYSIS



PART I: PATIENTS & PRESCRIPTIONS

METHODOLOGY

EXPLORATORY DATA ANALYSIS




PART 2: MEDICINES

METHODOLOGY

EXPLORATORY DATA ANALYSIS


Top 30 Prescribed Medicine



PART 2: MEDICINES

METHODOLOGY

EXPLORATORY DATA ANALYSIS



PART 2:
MEDICINES

METHODOLOGY

DATA PREPROCESSING & SPLITTING

STAGE 3



Prepare and ensure the quality of the data before used for model development

Missing Value

Resolve by skewness of data

Numeric Data

- if skewness ≈ 0 ; use mean
- if skewness < 0 or > 0 ; use median

Object Data

- use mode

Default Medicine

Map out top 50 ICD-10 with their popular prescribed medicine into json file.

- ICD name as key
- medicine & medicine count as value

Binning Patient Age

Age Range	Life Stage
0 - 12	Children
13 - 18	Adolescence
19 - 59	Adult
≥ 60	Senior Adult




Binning Patient BMI

Age Range	BMI Category
< 18.5	Underweight
18.5 - 24.9	Normal
≥ 25	Overweight

METHODOLOGY


AI MODEL DEVELOPMENT & EVALUATION

STAGE 4



MODEL 1: COLLABORATIVE FILTERING

"If a person A has the same opinion as a person B on an issue, A is more likely to have B's opinion on a different issue than that of a randomly chosen person."



MODEL 1: COLLABORATIVE FILTERING

SCORING METHOD


METHODOLOGY AI MODEL DEVELOPMENT & EVALUATION



Parameter	Range	Description	Score
ICD-10 Diagnosis	-	If all of the ICD-10 diagnoses match with the other patient, it will give a full score. If it only matches partially, ex. 2 out 5; hence it will be $\frac{2}{5}$.	4
Age Category	Difference: 1 - 2 [-0.25 Score] Difference: 2 - 4 [-0.50 Score] and so on [+ -0.25 Score]	For a different of more than one age category to the active patient, 0.25 will be subtracted from 1 (=100% similarity) for each additional age bin	1
BMI Category	-	If the BMI category matched, score is 1, else 0.	4
Gender	-	If gender matched score is 1, else 0.	1
Total Score			10

MODEL 1: COLLABORATIVE FILTERING

METHODOLOGY AI MODEL DEVELOPMENT & EVALUATION



MODEL 1: COLLABORATIVE FILTERING

INPUT DATA

Prescription	Prescription Pivot	Patients	Medicines	Medicine Prescriptions
Prescription ID	Patient ID	Patient ID	Medicine Brand	Frequency
Differential Diagnosis	Medicine Brand	Patient's Age Category	-	FrequencyDd
-	-	Patient's BMI Category	-	Timing
-	-	Gender	-	Duration
-	-	-	-	Amount

MODEL 1: COLLABORATIVE FILTERING

METHODOLOGY AI MODEL DEVELOPMENT & EVALUATION

OUTPUT DATA: MEDICINE LIST

```
{  
  "SURBEX T (KOMB Vit C 500 mg, Niasinamida 100 mg, Kalsium Pantotenat  
20 mg,Vit. B1 15 mg, Vit. B2 10 mg,Vit. B6 5 mg, Vit. B12 4 mg) TAB": 124,  
  "CETIRIZINE 10 MG Tablet": 78,  
  "PARACETAMOL 500 MG Tablet": 77,  
  "AMBROXOL 30 MG Tablet": 48,  
  "Acetylcystein 200mg Tablet": 33,  
  "SANADRYL DMP (Dekstrometorfan HBr 10 mg, Difenhidramin HCl 12,5 mg,  
Ammonium Cl 100 mg, Na Sitrat 50 mg, Mentol 1 mg/ 5ML) SYRUP, BOTOL  
120ML": 23,  
  "DEXAMETHASONE 0.5 MG Tablet": 20,  
  "PARACETAMOL 600 MG Tablet": 19,  
  "CETIRIZINE 5 MG/5 ML Botol": 18,  
  "ELKANA Syrup Botol 60ML": 18  
}
```

MODEL 1: COLLABORATIVE FILTERING


METHODOLOGY AI MODEL DEVELOPMENT & EVALUATION

OUTPUT DATA: MEDICINE DOSAGE

```
"SURBEX T (KOMB Vit C 500 mg, Niasinamida 100 mg, Kalsium Pantotenat 20  
mg,Vit. B1 15 mg, Vit. B2 10 mg,Vit. B6 5 mg, Vit. B12 4 mg) TAB": {  
  "frequency": 1.0,  
  "frequencyDd": 1.0,  
  "timing": "Setelah Makan",  
  "duration": 5.0,  
  "amount": 5.0  
}
```



MODEL 2: K-NEAREST NEIGHBOR

"Unknown Data point that is in a certain neighborhood will likely be similar with the most data point in that neighborhood"



MODEL 2: K-NEAREST NEIGHBOR

METHODOLOGY AI MODEL DEVELOPMENT & EVALUATION



MODEL 2: K-NEAREST NEIGHBOR

INPUT DATA

Prescription	Prescription Pivot	Patients	Medicines	Medicine Prescriptions
Prescription ID	Patient ID	Patient ID	Medicine Brand	Frequency
Differential Diagnosis	Medicine Brand	Patient's Age Category	-	FrequencyDd
-	-	Patient's BMI Category	-	Timing
-	-	Gender	-	Duration
-	-	-	-	Amount

MODEL 2: K-NEAREST NEIGHBOR

METHODOLOGY AI MODEL DEVELOPMENT & EVALUATION

OUTPUT DATA: MEDICINE LIST

```
{  
  "SURBEX T (KOMB Vit C 500 mg, Niasinamida 100 mg, Kalsium Pantotenat  
20 mg,Vit. B1 15 mg, Vit. B2 10 mg,Vit. B6 5 mg, Vit. B12 4 mg) TAB": 124,  
  "CETIRIZINE 10 MG Tablet": 78,  
  "PARACETAMOL 500 MG Tablet": 77,  
  "AMBROXOL 30 MG Tablet": 48,  
  "Acetylcystein 200mg Tablet": 33,  
  "SANADRYL DMP (Dekstrometorfan HBr 10 mg, Difenhidramin HCl 12,5 mg,  
Ammonium Cl 100 mg, Na Sitrat 50 mg, Mentol 1 mg/ 5ML) SYRUP, BOTOL  
120ML": 23,  
  "DEXAMETHASONE 0.5 MG Tablet": 20,  
  "PARACETAMOL 600 MG Tablet": 19,  
  "CETIRIZINE 5 MG/5 ML Botol": 18,  
  "ELKANA Syrup Botol 60ML": 18  
}
```

MODEL 2: K-NEAREST NEIGHBOR

METHODOLOGY AI MODEL DEVELOPMENT & EVALUATION



OUTPUT DATA: MEDICINE DOSAGE

```
"SURBEX T (KOMB Vit C 500 mg, Niasinamida 100 mg, Kalsium Pantotenat 20  
mg,Vit. B1 15 mg, Vit. B2 10 mg,Vit. B6 5 mg, Vit. B12 4 mg) TAB": {  
  "frequency": 1.0,  
  "frequencyDd": 1.0,  
  "timing": "Setelah Makan",  
  "duration": 5.0,  
  "amount": 5.0  
}
```

EVALUATION: PERFORMANCE METRICS

Recall@k:


The proportion of relevant medicines found in the recommendation medicines set

Precision@k:

Represents the proportion of recommended medicines (in the top k set; k=10) that are relevant

F1 Score:

Harmonic mean of both recall and precision



EVALUATION: K-FOLD CROSS VALIDATION



K Folds = 10
COLLABORATIVE FILTERING

Metric	Value
Recall@k	77.04
Precision@k	24.19
F1 Score	36.67
Runtime	2.14 s/item

K Folds = 10
K-NEAREST NEIGHBOR

Metric	Value
Recall@k	63.51
Precision@k	15.11
F1 Score	23.69
Runtime	1.06 s/item

MODEL REFINEMENT: BEST SCENARIOS



METHODOLOGY
MODEL REFINEMENT



COLLABORATIVE FILTERING

Threshold = 2.01

Metric	Value
Recall@k	66.30
Precision@k	23.69
F1 Score	33.78
Runtime	1.48 s/item

KNN

K ICD = 151 , K Patient = 51

Metric	Value
Recall@k	63.91
Precision@k	15.20
F1 Score	23.84
Runtime	1.07 s/item


KNN

K Patient = 150, K ICD = 40

Threshold: Similarity = 0.04 , Med = 0.4

Metric	Value
Recall@k	62.20
Precision@k	14.73
F1 Score	23.12
Runtime	1.34 s/item

MODEL REFINEMENT: BEST SCENARIOS



METHODOLOGY
MODEL REFINEMENT



COLLABORATIVE FILTERING

Threshold = 2.01

Metric	Value
Recall@k	66.30
Precision@k	23.69
F1 Score	33.78
Runtime	1.48 s/item

KNN

K ICD = 151 , K Patient = 51

Metric	Value
Recall@k	63.91
Precision@k	15.20
F1 Score	23.84
Runtime	1.07 s/item

KNN

K Patient = 150, K ICD = 40

Threshold: Similarity = 0.04 , Med = 0.4

Metric	Value
Recall@k	62.20
Precision@k	14.73
F1 Score	23.12
Runtime	1.34 s/item

MODEL DEPLOYMENT



User Input:

Patient Information:

- Patient ID
- Gender
- Age
- Height (cm)
- Weight (kg)

Differential Diagnosis

- Differential Diagnosis 1- 5

Welcome to e-Recommendation Prescription!

#Patient Information

We need some information to give the well-suited recommendations

Patient ID
c4d4a830-3277-11ea-a864-df576d4da157

Select Gender
Female

Slide your Age
1 35 150

Height (cm)
164 - +

Weight (kg)
58 - +

MODEL DEPLOYMENT



User Input:

Patient Information:

- Patient ID
- Gender
- Age
- Height (cm)
- Weight (kg)

Differential Diagnosis

- Differential Diagnosis 1- 5

Weight (kg)

58 - +

#Differential Diagnosis

Select Differential Diagnosis 1

R50.9 Fever, unspecified

Select Differential Diagnosis 2

K59.1 Functional diarrhoea

Select Differential Diagnosis 3

J00 Acute nasopharyngitis [common cold]

Select Differential Diagnosis 4

NaN

Select Differential Diagnosis 5

NaN

Give Recommendation

MEDICINE LIST

Give Recommendation

	0
patientId	c4d4a830-3277-11ea-a864-df576d4da157
patientGender	f
age	35
heights	164
weight	58
differentialDiagnosis	["R50.9 Fever, unspecified","K59.1 Functional diarrhoea","J00 Ac...
patientAgeCategory	Adult
patientBMICategory	Normal

▼ {

"PARACETAMOL 500 MG Tablet" : 4

"CETIRIZINE 5 MG/5 ML Botol" : 3

"PARACETAMOL 100 MG/ML Botol" : 3

"

SURBEX T (KOMB Vit C 500 mg, Niasinamida 100 mg, Kalsium Pantotenat 20 mg,Vit. B1 15 mg, Vit. B2 10 mg,Vit. B6 5 mg, Vit. B12 4 mg) TAB

"

: 2

"PARACETAMOL 600 MG Tablet" : 2

"PARACETAMOL 120 MG/5 ML Botol" : 1

"NEW DIATAB TABLET (MOLAGIT)" : 1

"DOMPERIDONE 10 MG Tablet" : 1

"CETIRIZINE 10 MG Tablet" : 1

"LOPERAMIDE 2 MG Tablet" : 1



METHODOLOGY
MODEL DEPLOYMENT

OUTPUT

MEDICINE DOSAGE

```

{
  "PARACETAMOL 500 MG Tablet" : {
    "frequency" : 3
    "frequencyDd" : 1
    "timing" : "Setelah Makan"
    "duration" : 3
    "amount" : 15
  }
  "CETIRIZINE 5 MG/5 ML Botol" : {
    "frequency" : 1
    "frequencyDd" : 5
    "timing" : "Setelah Makan"
    "duration" : 5
    "amount" : 1
  }
  "PARACETAMOL 100 MG/ML Botol" : {
    "frequency" : 3
    "frequencyDd" : 1
    "timing" : "Setelah Makan"
    "duration" : 5
    "amount" : 1
  }
  "SURBEX T (KOMB Vit C 500 mg, Niasinamida 100 mg, Kalsium Pantotenat 20 mg,Vit. B1 15 mg, Vit. B2 10 mg,Vit. B6 5 mg, Vit. B12 4 mg) TAB" :
  {

```



METHODOLOGY
MODEL DEPLOYMENT

OUTPUT



CONCLUSION

Best model developed is collaborative filtering with threshold 2.01. This model produce recall@k of 66.30, precision@k 23.69, f1 score 33.78, and runtime 1.48 s/item.





THANK YOU!

dave.peerera@student.sgu.ac.id

